

Tracheostomy tube with adjustable flange



- adjustable flange facilitate tube length adjustment
- made of thermosensitive PVC, siliconised
- with high volume low pressure cuff or without cuff
- X-ray line
- soft, smooth, transparent, adjustable flange
- obturator
- 2 neck tapes
- marked pilot balloon
- graduation mark
- latex free, phthalate free
- sterile, single use



size I.D. [mm]	O.D. [mm]	length L _{max} [mm]	reference number	
			with cuff	without cuff
5,0	6,9	58	35-5010	36-5010
5,5	7,5	62	35-5510	36-5510
6,0	8,2	72	35-6010	36-6010
6,5	8,9	76	35-6510	36-6510
7,0	9,7	84	35-7010	36-7010
7,5	10,3	92	35-7510	36-7510
8,0	11,0	97	35-8010	36-8010
8,5	11,7	102	35-8510	36-8510
9,0	12,3	104	35-9010	36-9010
9,5	13,0	109	35-9510	36-9510
10,0	13,7	110	35-1010	36-1010
11,0	14,8	110	35-1110	36-1110

Tracheostomy tube with adjustable flange, long



- extended proximal and distal length
- adjustable flange facilitate tube length adjustment
- made of thermosensitive PVC, siliconised
- with high volume low pressure cuff or without cuff
- X-ray line
- soft, smooth, transparent, adjustable flange
- obturator
- 2 neck tapes
- marked pilot balloon
- graduation mark
- latex free, phthalate free
- sterile, single use



size I.D. [mm]	O.D. [mm]	length L _{max} [mm]	reference number	
			with cuff	without cuff
7,0	9,7	100	35-7017	36-7017
8,0	11,0	116	35-8017	36-8017
9,0	12,3	122	35-9017	36-9017
10,0	13,7	128	35-1017	36-1017