

Microlaryngeal tracheal tube



- oral – nasal
- elongate length of the tube
- Murphy type
- made of thermosensitive PVC, siliconised
- transparent
- extra high volume low pressure cuff with thin walls
- X-ray line
- clear tube markings
- marked pilot balloon
- latex free, phthalate free
- sterile, single use



size I.D. [mm]	O.D. [mm]	Ø cuff [mm]	length [mm]	reference number	
				non-reinforced	reinforced
4,0	5,4	20,0	368	05-4010	05-4012
5,0	6,9	25,0	368	05-5010	05-5012
6,0	8,2	28,0	368	05-6010	05-6012