

## SUMINI neonate tracheostomy tube with cuff

New



- made of highest quality medical material, flexible and adapting to patient anatomy
- high volume low pressure cuff
- soft, flexible and special shaped flange facilitate firm fastening and minimizes skin irritation
- perfectly matched obturator
- neck tape for fastening
- visible in X-rays
- rounded, atraumatic tip
- latex free, phthalate free
- sterile, single use

size I.D. [mm]	O.D. [mm]	length L [mm]	Ø cuff [mm]	reference number
2,5	4,0	30	8,0	2N-2500
3,0	4,5	31	9,0	2N-3000
3,5	5,1	33	10,0	2N-3500
4,0	5,8	35	12,0	2N-4000

## SUMINI neonate tracheostomy tube without cuff

New



- made of highest quality medical material, flexible and adapting to patient anatomy
- soft, flexible and special shaped flange facilitate firm fastening and minimizes skin irritation
- perfectly matched obturator
- neck tape for fastening
- visible in X-rays
- rounded, atraumatic tip
- latex free, phthalate free
- sterile, single use

size I.D. [mm]	O.D. [mm]	length L [mm]	reference number
2,5	4,0	30	3N-2500
3,0	4,5	31	3N-3000
3,5	5,1	33	3N-3500
4,0	5,8	35	3N-4000