


# SURGICRYL<sup>®</sup>

910


POLYGLACTINE



<b>Type</b>	Braided and coated multifilament
<b>Composition</b>	Polyglactine 910, a copolymer made of 90% glycolide and 10%L-lactide
<b>Coating</b>	Poly-glycolide-co-L-lactide and calcium stearate (<1%)
<b>Colour</b>	Violet
<b>Tissue reaction</b>	Minimal
<b>Absorption</b>	The hydrolytic action by which the material is broken down results in total absorption between 56 and 70 days. Between 40% and 50% of tensile strength remain after 21 days.
<b>Presentation</b>	Box with 12 sutures



*For quality  
&  
safety*



[www.sutures.be](http://www.sutures.be)

# SURGICRYL®

## 910

### POLYGLACTINE

### CHARACTERISTICS

- High tensile strength
- Synthetic
- Breaks down by hydrolysis
- Braided structure
  
- Coated multifilament
  
- Hermetically sealed packing

### BENEFITS

- Excellent tissue support
- Minimal tissue reaction
- Predictable absorption rate
- Excellent handling properties
  - High flexibility
  - High tensile strength
  - Excellent knot security
- Soft passage through the tissues with minimal sawing, tissue drag and trauma
- Good knot tie down / adjustability
- Reduced capillary action
  
- Guaranteed seal and product sterility

**INDICATIONS** General soft tissue closing and/or ligation especially in


- General surgery
- Skin closure
- Gastrointestinal surgery
- Plastic surgery
- Ophthalmic surgery
- Gynaecology
- Obstetrics
- Orthopaedics
- Urology

**CONTRAINDICATIONS** Cardiovascular surgery  
Neurosurgery


Being absorbable not to use when an extended wound support is required.

Caution in patients with delayed wound healing (e.g. compromised immune system).

### TENSILE STRENGTH



### ABSORPTION



Without needle



Metric	USP	Length		Order number
--------	-----	--------	--	--------------

1.5	4/0	150 cm	<i>violet</i>	15150015
2	3/0	150 cm	<i>violet</i>	15200015
3	2/0	150 cm	<i>violet</i>	15300015
3.5	0	150 cm	<i>violet</i>	15350015
4	1	150 cm	<i>violet</i>	15400015
5	2	150 cm	<i>violet</i>	15500015
3	2/0	3 x 45 cm	<i>violet</i>	15300045-3
3.5	0	3 x 45 cm	<i>violet</i>	15350045-3
1.5	4/0	6 x 45 cm	<i>violet</i>	15150045-6
2	3/0	6 x 45 cm	<i>violet</i>	15200045-6
3	2/0	6 x 45 cm	<i>violet</i>	15300045-6
3.5	0	6 x 45 cm	<i>violet</i>	15350045-6
1.5	4/0	12 x 45 cm	<i>violet</i>	15150045
2	3/0	12 x 45 cm	<i>violet</i>	15200045
3	2/0	12 x 45 cm	<i>violet</i>	15300045
3.5	0	12 x 45 cm	<i>violet</i>	15350045
4	1	12 x 45 cm	<i>violet</i>	15400045
5	2	12 x 45 cm	<i>violet</i>	15500045
2	3/0	2 x 70 cm	<i>violet</i>	15200070-2
3	2/0	2 x 70 cm	<i>violet</i>	15300070-2
3.5	0	2 x 70 cm	<i>violet</i>	15350070-2
4	1	2 x 70 cm	<i>violet</i>	15400070-2
5	2	2 x 70 cm	<i>violet</i>	15500070-2

Without needle



Ligapack

















Metric	USP	Length		Order number
--------	-----	--------	--	--------------




2	3/0	250 cm	<i>violet</i>	15200025
3	2/0	250 cm	<i>violet</i>	15300025
3.5	0	250 cm	<i>violet</i>	15350025
4	1	250 cm	<i>violet</i>	15400025
5	2	250 cm	<i>violet</i>	15500025




Round bodied taper point 




1/2 Circle 






		Metric	USP	Length	Color	Order number
HR 13 mm		0.7	6/0	75 cm	<i>violet</i>	15070113
		1	5/0	75 cm	<i>violet</i>	15100113
		1.5	4/0	75 cm	<i>violet</i>	15150113
HR 17 mm		1	5/0	75 cm	<i>violet</i>	15100117
		1.5	4/0	75 cm	<i>violet</i>	15150117
		2	3/0	75 cm	<i>violet</i>	15200117
HR 22 mm		3	2/0	75 cm	<i>violet</i>	15300117
		1.5	4/0	75 cm	<i>violet</i>	15150122
		2	3/0	75 cm	<i>violet</i>	15200122
HR 26 mm		3	2/0	75 cm	<i>violet</i>	15300122
		4	1	75 cm	<i>violet</i>	15400122
		1.5	4/0	75 cm	<i>violet</i>	15150126
HR 30 mm		2	3/0	75 cm	<i>violet</i>	15200126
		3	2/0	75 cm	<i>violet</i>	15300126
		3.5	0	75 cm	<i>violet</i>	15350126
		4	1	75 cm	<i>violet</i>	15400126
		5	2	75 cm	<i>violet</i>	15500126
HR 36 mm		2	3/0	75 cm	<i>violet</i>	15200130
		3	2/0	75 cm	<i>violet</i>	15300130
		3.5	0	90 cm	<i>violet</i>	15350130
		4	1	90 cm	<i>violet</i>	15400130
		5	2	90 cm	<i>violet</i>	15500130
HR 40 mm		2	3/0	75 cm	<i>violet</i>	15200136
		3	2/0	75 cm	<i>violet</i>	15300136
		3.5	0	75 cm	<i>violet</i>	15350136
		3.5	0	90 cm	<i>violet</i>	15350140-75
		4	1	75 cm	<i>violet</i>	15400136
		4	1	90 cm	<i>violet</i>	15400140-75
		5	2	90 cm	<i>violet</i>	15500136
HR 45 mm		5	2	75 cm	<i>violet</i>	15500140-75
		5	2	90 cm	<i>violet</i>	15500140
		2	3/0	75 cm	<i>violet</i>	15200140
		3	2/0	75 cm	<i>violet</i>	15300140
		3.5	0	75 cm	<i>violet</i>	15350140-75












Round bodied taper point 		Metric	USP	Length		Order number
HR 48 mm 		3	2/0	90 cm	<i>violet</i>	15300148
		3.5	0	90 cm	<i>violet</i>	15350148
		4	1	90 cm	<i>violet</i>	15400148
		5	2	90 cm	<i>violet</i>	15500148
HR 60 mm 		3.5	0	75 cm	<i>violet</i>	15350160
		4	1	75 cm	<i>violet</i>	15400160
1/2 Circle  HR 65 mm 		3	2/0	90 cm	<i>violet</i>	15300165
		3.5	0	90 cm	<i>violet</i>	15350165
HR 80 mm 		4	1	100 cm	<i>violet</i>	15400180

Round bodied taper point double armed 		Metric	USP	Length		Order number
1/2 Circle  2 x HR 26 mm 		2	3/0	75 cm	<i>violet</i>	15200126-2

Round bodied taper point loop 		Metric	USP	Length		Order number
1/2 Circle  HR 45 mm 		3.5	0	150 cm	<i>violet</i>	15350945

Round bodied taper point 		Metric	USP	Length		Order number
<b>5/8 Circle</b> 	FR 27 mm 	3	2/0	75 cm	<i>violet</i>	15300327
		3.5	0	75 cm	<i>violet</i>	15350327

Round bodied taper point 		Metric	USP	Length		Order number
<b>3/8 Circle</b> 	DR 17 mm 	1.5	4/0	45 cm	<i>violet</i>	15150217
	DR 22 mm 	2	3/0	75 cm	<i>violet</i>	15200222
	DR 30 mm 	4	1	75 cm	<i>violet</i>	15400230

Reverse cutting 		Metric	USP	Length		Order number
<b>1/2 Circle</b> 	HS 12 mm 	0.7	6/0	45 cm	<i>violet</i>	15071412
	HS 17 mm 	1	5/0	75 cm	<i>violet</i>	15101417
		1.5	4/0	75 cm	<i>violet</i>	15151417
	HS 22 mm 	2	3/0	75 cm	<i>violet</i>	15201422
	HS 23 mm 	5	2	90 cm	<i>violet</i>	15501423
	HS 26 mm 	2	3/0	75 cm	<i>violet</i>	15201426
		3	2/0	75 cm	<i>violet</i>	15301426
		3.5	0	75 cm	<i>violet</i>	15351426
		4	1	75 cm	<i>violet</i>	15401426
	HS 30 mm 	2	3/0	75 cm	<i>violet</i>	15201430
	HS 36 mm 	3	2/0	75 cm	<i>violet</i>	15301436
		3.5	0	90 cm	<i>violet</i>	15351436
		4	1	90 cm	<i>violet</i>	15401436
		5	2	90 cm	<i>violet</i>	15501436
	HS 40 mm 	2	3/0	90 cm	<i>violet</i>	15201440
3		2/0	90 cm	<i>violet</i>	15301440	
3.5		0	90 cm	<i>violet</i>	15351440	
4		1	90 cm	<i>violet</i>	15401440	
HS 48 mm 	5	2	90 cm	<i>violet</i>	15501440	
	3.5	0	90 cm	<i>violet</i>	15351448	
	4	1	90 cm	<i>violet</i>	15401448	
	5	2	90 cm	<i>violet</i>	15501448	



Reverse cutting
















3/8 Circle











		Metric	USP	Length		Order number
DS 12 mm		0.7	6/0	45 cm	<i>violet</i>	15071512
		1	5/0	45 cm	<i>violet</i>	15101512
		1.5	4/0	45 cm	<i>violet</i>	15151512
DS 16 mm		0.7	6/0	45 cm	<i>violet</i>	15071516
		1	5/0	75 cm	<i>violet</i>	15101516
		1.5	4/0	75 cm	<i>violet</i>	15151516
		2	3/0	45 cm	<i>violet</i>	15201516-45
DS 16 mm Easy Pass		2	3/0	50 cm	<i>undyed</i>	15201516UP
DS 16 mm black		0.7	6/0	45 cm	<i>violet</i>	15071516VB
DS 19 mm		1	5/0	75 cm	<i>violet</i>	15101519
		1.5	4/0	75 cm	<i>violet</i>	15151519
		2	3/0	75 cm	<i>violet</i>	15201519
DS 19 mm Easy Pass		1.5	4/0	50 cm	<i>undyed</i>	15151519UP
DS 19 mm black		1	5/0	75 cm	<i>violet</i>	15101519VB
		1.5	4/0	75 cm	<i>violet</i>	15151519VB
DS 24 mm		1.5	4/0	75 cm	<i>violet</i>	15151524
		2	3/0	75 cm	<i>violet</i>	15201524
		3	2/0	75 cm	<i>violet</i>	15301524
		3.5	0	75 cm	<i>violet</i>	15351524
		DS 24 mm Easy Pass		2	3/0	50 cm
DS 26 mm		2	3/0	75 cm	<i>violet</i>	15201526
		3	2/0	75 cm	<i>violet</i>	15301526
DS 26 mm Easy Pass		2	3/0	75 cm	<i>undyed</i>	15201526UP
DS 30 mm		2	3/0	75 cm	<i>violet</i>	15201530
		3	2/0	75 cm	<i>violet</i>	15301530
		3.5	0	75 cm	<i>violet</i>	15351530
DS 36 mm		3	2/0	75 cm	<i>violet</i>	15301536
		3.5	0	75 cm	<i>violet</i>	15351536
		4	1	75 cm	<i>violet</i>	15401536
DS 40 mm		3.5	0	75 cm	<i>violet</i>	15351540
		4	1	75 cm	<i>violet</i>	15401540
DS 45 mm		3	2/0	75 cm	<i>violet</i>	15301545


Reverse cutting 	Metric	USP	Length	Order number	
<b>3/8 circle</b> 	DS 80 mm	7	5	150 cm	<i>violet</i> 15701580


Round bodied taper cutting point 	Metric	USP	Length	Order number	
<b>1/2 Circle</b> 	HRT 13 mm 	1.5	4/0	75 cm	<i>violet</i> 15150513
	HRT 17 mm 	1.5	4/0	75 cm	<i>violet</i> 15150517
		2	3/0	75 cm	<i>violet</i> 15200517
	HRT 26 mm 	1.5	4/0	75 cm	<i>violet</i> 15150526
		2	3/0	75 cm	<i>violet</i> 15200526
	HRT 30 mm 	3	2/0	75 cm	<i>violet</i> 15300526
		3	2/0	75 cm	<i>violet</i> 15300530
	HRT 36 mm 	3.5	0	75 cm	<i>violet</i> 15350530
		4	1	75 cm	<i>violet</i> 15400530
	HRT 40 mm 	3	2/0	90 cm	<i>violet</i> 15300536
		3.5	0	90 cm	<i>violet</i> 15350536
	HRT 48 mm 	4	1	90 cm	<i>violet</i> 15400536
4		1	90 cm	<i>violet</i> 15400540	
	5	2	90 cm	<i>violet</i> 15500540	
	8	6	150 cm	<i>violet</i> 15800548	


Round bodied taper cutting point 	Metric	USP	Length	Order number	
<b>3/8 Circle</b> 	DRT 13 mm 	1.5	4/0	75 cm	<i>violet</i> 15150413
	DRT 17 mm 	1.5	4/0	75 cm	<i>violet</i> 15150417
		2	3/0	75 cm	<i>violet</i> 15200417



Round bodied taper cutting point 		Metric	USP	Length	Order number
<b>3/8 Circle</b>  <b>1/2 Circle</b> 	DRT 13 mm	1	5/0	75 cm	<i>violet</i> 15100530-0413
	HRT 30 mm				
Cutting 		Metric	USP	Length	Order number
<b>Straight</b>	GS 60 mm	2	3/0	75 cm	<i>violet</i> 15202160
		3	2/0	75 cm	<i>violet</i> 15302160
Blunt 		Metric	USP	Length	Order number
<b>1/2 Circle</b> 	HRB 30 mm	3	2/0	75 cm	<i>violet</i> 15301830
	HRB 36 mm	3.5	0	90 cm	<i>violet</i> 15351836
		4	1	90 cm	<i>violet</i> 15401836
HRB 45 mm	5	2	90 cm	<i>violet</i> 15501845	
Blunt 		Metric	USP	Length	Order number
<b>3/8 Circle</b> 	DRB 63 mm	4	1	100 cm	<i>violet</i> 15401963
	DRB 85 mm	5	2	75 cm	<i>violet</i> 15501985
Micro needle taper point 		Metric	USP	Length	Order number
<b>3/8 Circle</b> 	MDR 6,4 mm 150µm	0.4	8/0	20 cm	<i>violet</i> M910087

Side cutting spatula double armed		Metric	USP	Length	Order number	
<b>1/2 Circle</b> 	2 x HSP 8 mm 200 µm	0.4	8/0	45 cm	<i>violet</i>	O910081
	2 x HSP 8 mm 200 µm	0.5	7/0	45 cm	<i>violet</i>	O91007
	300 µm	0.7	6/0	45 cm	<i>violet</i>	O910061
	2 x HSP 8 mm 300 µm	0.7	6/0	45 cm	<i>undyed</i>	O910061U

Side cutting spatula double armed		Metric	USP	Length	Order number	
<b>3/8 Circle</b> 	2 x DSP 6 mm 200 µm	0.4	8/0	45 cm	<i>violet</i>	O91008
	220 µm	0.5	7/0	45 cm	<i>violet</i>	O910072
	2 x DSP 11 mm 350 µm	1	5/0	45 cm	<i>violet</i>	O910051

Side cutting spatula double armed		Metric	USP	Length	Order number	
<b>1/4 circle</b> 	2 x VSP 8 mm 300 µm	0.7	6/0	45 cm	<i>violet</i>	O91006
		0.7	6/0	20 cm	<i>undyed</i>	O91006U-20
		0.7	6/0	45 cm	<i>undyed</i>	O91006U
	2 x VSP 8 mm 350 µm	1	5/0	45 cm	<i>violet</i>	O91005
		1	5/0	45 cm	<i>undyed</i>	O91005U