


Closed suction system 72h



- mechanical push block system, sealing off access to the patient
- double swivel elbow connector
- MDI port for quick and effective drug delivery
- one-way valve for catheter rinsing
- transparent and durable sleeve securing suction catheter
- catheter with clearly visible graduation and rounded tip to minimize the risk of trachea mucosa damaging
- double black markings of catheter withdrawal
- colour coded sizes of suction catheter
- suction control button with built-in locking system to prevent accidental pressing
- sterile, phthalate free, latex free

Set includes

- closed suction catheter
- extendible tube
- disconnecting wedge
- day sticker

size [CH/FG]	colour	tracheal		tracheostomy	
		length [mm]	reference	length [mm]	reference
10	Black	540	73-1071	360	74-1071
12	White	540	73-1271	360	74-1271
14	Green	540	73-1471	360	74-1471
16	Orange	540	73-1671	360	74-1671
18	Red	540	73-1871	360	74-1871

Closed suction system 24h



- double swivel elbow connector
- MDI port for quick and effective drug delivery
- one-way valve for catheter rinsing
- catheter with clearly visible graduation and rounded tip to minimize the risk of trachea mucosa damaging
- transparent and durable sleeve securing suction catheter
- double black markings of catheter withdrawal
- colour coded sizes of suction catheter
- suction control button with built-in locking system to prevent accidental pressing
- sterile, phthalate free, latex free

Set includes

- closed suction catheter
- extendible tube
- disconnecting wedge
- day sticker

size [CH/FG]	colour	tracheal	
		length [mm]	reference
10	Black	540	73-1023
12	White	540	73-1223
14	Green	540	73-1423
16	Orange	540	73-1623
18	Red	540	73-1823