

## Double lumen bronchial tube - right sided



- made of thermosensitive PVC, siliconised
- two delicate high volume low pressure cuffs
- smooth walls preventing gathering of secretion
- atraumatic, rounded tip and smooth cuff welding
- X-ray line along the tube, additional X-ray marks assisting in accurate positioning of the tube
- stylet covered with medical plastic
- marked pilot balloons
- graduation mark
- two 15 mm connectors conforming to ISO 5356-1
- latex free, phthalate free
- sterile, single use

size [CH/FG]	Ø <sub>B</sub> cuff [mm]	Ø <sub>T</sub> cuff [mm]	reference number	
			-	with carina hook
<b>24</b>	12,0	14,0	22-2410	-
<b>26</b>	15,0	20,0	22-2610	-
<b>28</b>	16,0	20,0	22-2810	-
<b>32</b>	16,0	21,0	22-3210	22-3241
<b>35</b>	18,0	23,0	22-3510	22-3541
<b>37</b>	18,0	25,0	22-3710	22-3741
<b>39</b>	20,0	30,0	22-3910	22-3941
<b>41</b>	22,0	32,0	22-4110	22-4141

## Double lumen bronchial tube – right sided, set



Set contains:

- bronchial tube, right sided
- connectors set, with or without clamps
- 4 control valve suction catheters

size [CH/FG]	Ø <sub>B</sub> cuff [mm]	Ø <sub>T</sub> cuff [mm]	reference number		
			bronchial tube		
			connectors without clamps	connectors with clamps	carina hook, connectors with clamps
<b>24</b>	12,0	14,0	22-2420	22-2422	-
<b>26</b>	15,0	20,0	22-2620	22-2622	-
<b>28</b>	16,0	20,0	22-2820	22-2822	-
<b>32</b>	16,0	21,0	22-3220	22-3222	22-3245
<b>35</b>	18,0	23,0	22-3520	22-3522	22-3545
<b>37</b>	18,0	25,0	22-3720	22-3722	22-3745
<b>39</b>	20,0	30,0	22-3920	22-3922	22-3945
<b>41</b>	22,0	32,0	22-4120	22-4122	22-4145