


KAN tracheostomy tube without cuff, fenestrated

Set I



Set components:


- KAN tube without cuff, fenestrated (1)
- obturator (2)
- 1 inner cannula with 15 mm connector fenestrated (3)
- 1 inner cannula with 15 mm connector, non-fenestrated (4)
- 1 inner cannula with collar, fenestrated (5)
- soft tracheostomy tube holder (6)
- 15 mm connector (7)
- speaking valve (8)
- connecting plug (9)
- decannulation cap (10)
- protective plug (11)

size I.D. [mm]	colour	tube O.D. [mm]	cannula I.D. [mm]	length L [mm]	tube I.D. ₁ [mm]	bending angle Θ [°]	ref.
5,0	Yellow	9,1	5,0	62	6,5	115	47-5010
6,0	Orange	10,2	6,0	70	7,5	110	47-6010
7,0	Green	11,2	7,0	72	8,5	105	47-7010
8,0	White	12,2	8,0	80	9,5	100	47-8010
9,0	Blue	13,1	9,0	82	10,5	95	47-9010
10,0	Yellow	14,3	10,0	84	11,5	90	47-1010

- KAN tube made of thermosensitive PVC
- perfectly matched transparent tube and inner cannulas
- concentric fenestrated holes
- soft flanges with joint hitch
- bayonet joint
- colour coded
- facilitates cleaning and disinfection of cannulas
- obturator
- phthalate free, latex free
- sterile

KAN tracheostomy tube without cuff, fenestrated

Set II



Set components:

- KAN tube without cuff, fenestrated (1)
- obturator (2)
- 2 inner cannulas with collar, fenestrated (3)
- soft tracheostomy tube holder (4)
- speaking valve (5)

size I.D. [mm]	colour	tube O.D. [mm]	cannula I.D. [mm]	length L [mm]	tube I.D. ₁ [mm]	bending angle Θ [°]	ref.
5,0	Yellow	9,1	5,0	62	6,5	115	47-5011
6,0	Orange	10,2	6,0	70	7,5	110	47-6011
7,0	Green	11,2	7,0	72	8,5	105	47-7011
8,0	White	12,2	8,0	80	9,5	100	47-8011
9,0	Blue	13,1	9,0	82	10,5	95	47-9011
10,0	Yellow	14,3	10,0	84	11,5	90	47-1011