

Preformed tracheal tube, nasal, soft



- nasal
- Murphy type
- made of thermosensitive PVC, siliconised
- milky, non-transparent
- with high volume low pressure cuff or without cuff
- flexible, with thicker walls, non kinking
- curved towards the forehead
- designated for surgeries in head and neck area
- X-ray line
- clear tube markings
- marked pilot balloon
- latex free, phthalate free
- sterile, single use

size I.D. [mm]	O.D. [mm]	Ø cuff [mm]	reference number	
			with cuff	without cuff
3,0	4,5	8,0	27-3014	28-3014
3,5	5,1	10,0	27-3514	28-3514
4,0	6,0	10,0	27-4014	28-4014
4,5	6,7	12,0	27-4514	28-4514
5,0	7,3	14,0	27-5014	28-5014
5,5	8,0	17,0	27-5514	28-5514
6,0	8,8	20,0	27-6014	28-6014
6,5	9,4	20,0	27-6514	28-6514
7,0	10,2	23,0	27-7014	28-7014
7,5	10,8	25,0	27-7514	28-7514
8,0	11,6	26,0	27-8014	28-8014

Preformed tracheal tube, nasal



- nasal
- Murphy type
- made of thermosensitive PVC, siliconised
- transparent
- with high volume low pressure cuff or without cuff
- extra soft
- curved towards the forehead
- designated for surgeries in head area
- X-ray line
- clear tube markings
- marked pilot balloon
- latex free, phthalate free
- sterile, single use

size I.D. [mm]	O.D. [mm]	Ø cuff [mm]	reference number	
			with cuff	without cuff
3,0	4,2	8,0	27-3010	28-3010
3,5	4,8	10,0	27-3510	28-3510
4,0	5,4	10,0	27-4010	28-4010
4,5	6,0	12,0	27-4510	28-4510
5,0	6,9	14,0	27-5010	28-5010
5,5	7,5	17,0	27-5510	28-5510
6,0	8,2	20,0	27-6010	28-6010
6,5	8,8	20,0	27-6510	28-6510
7,0	9,7	23,0	27-7010	28-7010
7,5	10,2	25,0	27-7510	28-7510
8,0	11,0	26,0	27-8010	28-8010