



Preformed tracheal tube, oral



- oral
- Murphy type
- made of thermosensitive PVC, siliconised
- transparent
- with high volume low pressure cuff or without cuff
- curved towards the chin
- designated for surgeries in head area
- X-ray line
- clear tube markings
- marked pilot balloon
- latex free
- sterile, single use

size I.D. [mm]	O.D. [mm]	Ø cuff [mm]	reference number	
			with cuff	without cuff
3,0	4,2	8,0	25-3010	26-3010
3,5	4,8	10,0	25-3510	26-3510
4,0	5,4	10,0	25-4010	26-4010
4,5	6,0	12,0	25-4510	26-4510
5,0	6,9	14,0	25-5010	26-5010
5,5	7,5	20,0	25-5510	26-5510
6,0	8,2	20,0	25-6010	26-6010
6,5	8,8	20,0	25-6510	26-6510
7,0	9,7	23,0	25-7010	26-7010
7,5	10,2	25,0	25-7510	26-7510
8,0	11,0	26,0	25-8010	26-8010
8,5	11,6	27,0	25-8510	26-8510
9,0	12,3	29,0	25-9010	26-9010

Preformed tracheal tube, oral – forehead



- oral
- Murphy type
- made of thermosensitive PVC, siliconised
- transparent
- with high volume low pressure cuff or without cuff
- curved towards the forehead
- designated for surgeries in head area
- X-ray line
- clear tube markings
- marked pilot balloon
- latex free, phthalate free
- sterile, single use

size I.D. [mm]	O.D. [mm]	Ø cuff [mm]	reference number	
			with cuff	without cuff
3,0	4,2	-	-	26-3037
3,5	4,8	-	-	26-3537
4,0	5,4	-	-	26-4037
4,5	6,0	-	-	26-4537
5,0	6,9	14,0	25-5037	26-5037
5,5	7,5	20,0	25-5537	26-5537
6,0	8,2	20,0	25-6037	26-6037
6,5	8,8	20,0	25-6537	26-6537
7,0	9,7	23,0	25-7037	26-7037
7,5	10,3	25,0	25-7537	26-7537
8,0	11,0	26,0	25-8037	26-8037
8,5	11,7	27,0	25-8537	26-8537
9,0	12,3	29,0	25-9037	26-9037