

## SUMINI paediatric tracheostomy tube with cuff

New



- made of highest quality medical material, flexible and adapting to patient anatomy
- high volume low pressure cuff
- soft, flexible and special shaped flange facilitate firm fastening and minimizes skin irritation
- perfectly matched obturator
- neck tape for fastening
- visible in X-rays
- rounded, atraumatic tip
- latex free, phthalate free
- sterile, single use

size I.D. [mm]	O.D. [mm]	length L [mm]	Ø cuff [mm]	reference number
2,5	4,0	32	8,0	2P-2500
3,0	4,5	35	9,0	2P-3000
3,5	5,1	39	10,0	2P-3500
4,0	5,8	43	12,0	2P-4000
4,5	6,5	47	12,0	2P-4500
5,0	7,1	49	15,0	2P-5000
5,5	7,6	54	15,0	2P-5500
6,0	8,3	60	18,0	2P-6000

## SUMINI paediatric tracheostomy tube without cuff

New



- made of highest quality medical material, flexible and adapting to patient anatomy
- soft, flexible and special shaped flange facilitate firm fastening and minimizes skin irritation
- perfectly matched obturator
- neck tape for fastening
- visible in X-rays
- rounded, atraumatic tip
- latex free, phthalate free
- sterile, single use

size I.D. [mm]	O.D. [mm]	length L [mm]	reference number
2,5	4,0	32	3P-2500
3,0	4,5	35	3P-3000
3,5	5,1	39	3P-3500
4,0	5,8	43	3P-4000
4,5	6,5	47	3P-4500
5,0	7,1	49	3P-5000
5,5	7,6	54	3P-5500
6,0	8,3	60	3P-6000