

Silicone laryngeal mask, reusable



- made of medical grade silicone
- transparent tube and soft silicone cuff
- anatomic shape for easy, gentle insertion and optimal placement
- size, weight and cuff volume indicated on tube
- 15 mm connector
- reusable up to 40 times
- latex free, phthalate free
- sterile

size	patient weight [kg]	cuff volume [ml]	reference number
1,0	0-5	4	42-1000
1,5	5-10	7	42-1500
2,0	10-20	10	42-2000
2,5	20-30	14	42-2500
3,0	30-50	20	42-3000
4,0	50-70	30	42-4000
5,0	70-100	40	42-5000

Silicone laryngeal mask, single use



- made of medical grade silicone
- transparent tube and soft silicone cuff
- anatomic shape for easy, gentle insertion and optimal placement
- size, weight and cuff volume indicated on tube
- colour coded
- 15 mm connector
- latex free, phthalate free
- sterile, single use

size	colour	patient weight [kg]	cuff volume [ml]	reference number	
				standard	reinforced
1,0		0-5	4	42-1010	42-1012
1,5		5-10	7	42-1510	42-1512
2,0		10-20	10	42-2010	42-2012
2,5		20-30	14	42-2510	42-2512
3,0		30-50	20	42-3010	42-3012
4,0		50-70	30	42-4010	42-4012
5,0		70-100	40	42-5010	42-5012

Available in set with syringe and lubricating jelly under reference number with „J” at the end, for example „42-4010J.”

Laryngeal mask, single use



- made of medical grade PVC
- transparent tube and soft cuff
- anatomic shape for easy, gentle insertion and optimal placement
- size, weight and cuff volume indicated on tube
- colour coded
- 15 mm connector
- latex free, phthalate free
- sterile, single use

size	colour	patient weight [kg]	cuff volume [ml]	reference number
1,0		0-5	4	42-1050
1,5		5-10	7	42-1550
2,0		10-20	10	42-2050
2,5		20-30	14	42-2550
3,0		30-50	20	42-3050
4,0		50-70	30	42-4050
5,0		70-100	40	42-5050

Available in set with syringe and lubricating jelly under reference number with „J” at the end, for example „42-4050J.”