

Tracheostomy tube with suction lumen



- made of thermosensitive PVC, siliconised
- very usefull for long-term ventilation
- allows suctioning of secretions from above the cuff
- suction lumen ended with universal connector
- high volume low pressure cuff
- X-ray line
- soft, smooth, transparent flange - fixed or adjustable
- obturator
- 2 neck tapes
- marked pilot balloon
- latex free, phthalate free
- sterile, single use



size I.D. [mm]	O.D. [mm]	reference number		reference number	
		length L [mm]	fixed flange	length L _{max} [mm]	adjustable flange
5,0	6,9	48	39-5041	58	39-5011
5,5	7,5	50	39-5541	62	39-5511
6,0	8,2	59	39-6041	72	39-6011
6,5	8,9	62	39-6541	77	39-6511
7,0	9,7	68	39-7041	85	39-7011
7,5	10,3	74	39-7541	93	39-7511
8,0	11,0	82	39-8041	98	39-8011
8,5	11,7	84	39-8541	102	39-8511
9,0	12,3	93	39-9041	105	39-9011
9,5	13,0	97	39-9541	110	39-9511
10,0	13,7	98	39-1041	111	39-1011