

Tracheal intubation stylet

- metal alloy covered with medical plastic
- soft tip
- latex free, phthalate free
- sterile, single use

size O.D. [mm]	length [mm]	reference number	for tubes size I.D. [mm]
1,9	230	09-1923	2,0 ÷ 3,0
2,0	230	09-2023	2,5 ÷ 3,0
2,2	230	09-2223	2,5 ÷ 3,5
3,0	340	09-3034	3,5 ÷ 5,0
4,0	340	09-4034	5,0 ÷ 8,0
4,0	600	09-4060	5,0 ÷ 8,0
5,0	370	09-5037	6,0 ÷ 11,0
5,0	600	09-5060	6,0 ÷ 11,0

