

Tracheal tube with low pressure cuff



- oral – nasal
- Murphy type
- made of thermosensitive PVC
- transparent
- high volume low pressure cuff
- X-ray line
- clear tube markings
- marked pilot balloon
- transparent 15 mm connector
- wide black ring depth mark
- latex free
- sterile, single use

size I.D. [mm]	O.D. [mm]	Ø cuff [mm]	reference number	
			PVC	PVC siliconised
3,0	4,2	12,0	02-3000	02-3010
3,5	4,9	12,0	02-3500	02-3510
4,0	5,5	14,0	02-4000	02-4010
4,5	6,2	14,0	02-4500	02-4510
5,0	6,9	17,0	02-5000	02-5010
5,5	7,5	17,0	02-5500	02-5510
6,0	8,2	20,0	02-6000	02-6010
6,5	8,8	20,0	02-6500	02-6510
7,0	9,6	25,0	02-7000	02-7010
7,5	10,2	25,0	02-7500	02-7510
8,0	10,9	26,0	02-8000	02-8010
8,5	11,5	26,0	02-8500	02-8510
9,0	12,1	28,0	02-9000	02-9010
9,5	12,7	28,0	02-9500	02-9510
10,0	13,6	28,0	02-1000	02-1010

Tracheal tube with low pressure cuff



- improved tube flexibility
- oral – nasal
- Murphy type
- made of thermosensitive PVC, siliconised
- transparent
- high volume low pressure cuff
- X-ray line
- clear tube markings
- marked pilot balloon
- different types of depth markings, available
- latex free, phthalate free
- sterile, single use

size I.D. [mm]	O.D. [mm]	Ø cuff [mm]	reference number
2,0	3,0	7,0	02-2012
2,5	3,5	8,0	02-2512
3,0	4,2	8,0	02-3012
3,5	4,8	10,0	02-3512
4,0	5,4	10,0	02-4012
4,5	6,0	12,0	02-4512
5,0	6,9	14,0	02-5012
5,5	7,5	17,0	02-5512
6,0	8,2	20,0	02-6012
6,5	8,8	20,0	02-6512
7,0	9,7	23,0	02-7012
7,5	10,3	25,0	02-7512
8,0	11,0	26,0	02-8012
8,5	11,7	27,0	02-8512
9,0	12,3	29,0	02-9012
9,5	13,0	29,0	02-9512
10,0	13,7	32,0	02-1012