

Tracheal tube with suction lumen



- oral - nasal
- Murphy type
- made of thermosensitive PVC, siliconised
- transparent
- high volume low pressure cuff
- allows suctioning of secretions from above the cuff, minimizing the risk of infection
- suction lumen ended with universal connector
- X-ray line
- clear tube markings
- marked pilot balloon
- recommended for long – term ventilation
- latex free, phthalate free
- sterile, single use

size I.D. [mm]	O.D. [mm]	Ø cuff [mm]	reference number	
			without stylet	with stylet
5,0	7,6	20,0	06-5011	06-5014
5,5	8,3	23,0	06-5511	06-5514
6,0	9,0	23,0	06-6011	06-6014
6,5	9,7	23,0	06-6511	06-6514
7,0	10,4	26,0	06-7011	06-7014
7,5	11,0	27,0	06-7511	06-7514
8,0	11,7	30,0	06-8011	06-8014
8,5	12,2	32,0	06-8511	06-8514
9,0	12,8	36,0	06-9011	06-9014
9,5	13,4	36,0	06-9511	06-9514
10,0	14,2	36,0	06-1011	06-1014