

Tracheostomy tube with high volume low pressure cuff



- made of thermosensitive PVC, siliconised
- high volume low pressure cuff
- with or without fenestration
- X-ray line
- soft, smooth, transparent flange
- obturator
- 2 neck tapes
- marked pilot balloon
- latex free, phthalate free
- sterile, single use

size I.D. [mm]	O.D. [mm]	length L [mm]	reference number	
			standard	fenestrated
3,0	4,2	35	31-3010	-
3,5	4,8	37	31-3510	-
4,0	5,4	40	31-4010	-
4,5	6,0	46	31-4510	-
5,0	6,9	48	31-5010	31-5055
5,5	7,5	50	31-5510	31-5555
6,0	8,2	59	31-6010	31-6055
6,5	8,9	62	31-6510	31-6555
7,0	9,7	68	31-7010	31-7055
7,5	10,3	74	31-7510	31-7555
8,0	11,0	82	31-8010	31-8055
8,5	11,7	84	31-8510	31-8555
9,0	12,3	93	31-9010	31-9055
9,5	13,0	97	31-9510	31-9555
10,0	13,7	98	31-1010	31-1055

Tracheostomy tube with double cuff



- made of thermosensitive PVC, siliconised
- two low pressure cuffs
- X-ray line
- soft, smooth, transparent flange - fixed or adjustable
- obturator
- 2 neck tapes
- marked pilot balloon
- latex free, phthalate free
- sterile, single use

size I.D. [mm]	O.D. [mm]	reference number		reference number	
		length L [mm]	with fixed flange	length L _{max} [mm]	with adjustable flange
5,0	6,9	58	32-5010	68	32-5040
6,0	8,2	67	32-6010	80	32-6040
7,0	9,7	79	32-7010	99	32-7040
7,5	10,3	82	32-7510	106	32-7540
8,0	11,0	92	32-8010	113	32-8040
8,5	11,7	92	32-8510	116	32-8540
9,0	12,3	107	32-9010	118	32-9040
9,5	13,0	109	32-9510	120	32-9540
10,0	13,7	115	32-1010	126	32-1040